

**Saudi Mais Co., for Medical Products**  
**Work Instruction**

**Title**

**Production & Test Procedure for Guedel Air Way**

|        |   |             |
|--------|---|-------------|
| Doc. # | : | WI-8-2-2-27 |
| Rev. # | : | 7           |
| Date.# | : | 29.08.2020  |
| Page.# | : | 1 / 5       |

**1. Purpose:**

To stipulate the required tests for Guedel air way to ensure good quality products before releasing to the market.

**2. Scope:**

Medical Device – Guedel Air Way.

**3. Definition:**

**Guedel Air Way** : A medical device used for used to give access to or pharyngeal pathway during or following general anesthesia.

**4. Responsibilities:**

- 4.1 Quality Controller is responsible for sampling & testing.
- 4.2 Quality control supervisor is responsible for checking the obtained results.
- 4.3 Machine operators and assembly peoples are responsible to produced good quantity and good quality products.

**5. Procedure:**

**5.1 Receiving Raw material:**

- Ensure the receiving slip from supplier is correct material according to specification.
- Inspect the material bag condition, must be no damaged at all bags and without any irregularities and dust particles.
- Take samples & inspect in order to AQL procedure, product description, specification with apply the AQL sampling plans for normal inspection.
- Record into OBM inspection report.

**5.1.1 Inspection & test procedure for Raw Material:**

| Work Description     | Requirements            | Method of Check |
|----------------------|-------------------------|-----------------|
| Clarity / Appearance | Clear / Good Appearance | Visual          |
| Condition of Bag     | Good / No damage bags   | Visual          |

Prepared by: QA Coordinator  
M.Azeem

Approved by: Q.C.M  
Mr. Javeed Gafoor

**Controlled Copy**

**Li Mais Co., for Medical Products**  
**Work Instruction**

**Production & Test Procedure for Guedel Air Way**

|        |   |             |
|--------|---|-------------|
| Doc. # | : | WI-8-2-2-27 |
| Rev. # | : | 7           |
| Date.# | : | 29.08.2020  |
| Page.# | : | 2 / 5       |

| Material Grades / Description | As per P.O                      | Visual / Physical |
|-------------------------------|---------------------------------|-------------------|
| Appearance of granules        | Good / Transparent / Colored    | Visual            |
| Granules formation            | No lobes formation              | Visual            |
| Lot number printed status     | Clear printed                   | Visual            |
| Arrangement of materials bags | Should be specified clear label | Visual            |

**5.1 Injection Machine Procedure:**

- 5.1.1 Production controller is responsible to issue required raw material for airway body and bite block and issue production order after receiving raw material.
- 5.1.2 Mold maintenance technician is responsible to prepare required mold for Guedel air way parts.
- 5.1.3 Machine maintenance technicians are responsible to mount the mold on specified machine and set approved parameters.
- 5.1.4 Machine operator is responsible to make sure that material is same as per production order and feed the material into mutant feeding machine or feeding through by manual.
- 5.1.5 Machine operator is responsible to give initial sample to quality controller for evaluation, quality controller is responsible to evaluate the products & make sure that the product is meeting the criteria of specification with the following quality inspection points.
- 5.1.5.1 Airway body**  
Quality controller is responsible to make sure that airway body accurate in dimension with smooth surface, no flashing, no black spots, no burn marks, flange must be round, without sharp edges, without dust particles and spots and record the results in process control form.
- 5.1.5.2 Bite Block**  
Quality controller is responsible to make sure that bite block must be accurate in dimension, color coded, smooth surface, no flashing, no black spot, no sharp edges and free from irregularities defects and record the results in process control form.
- 5.1.6 Machine Maintenance technicians are responsible to adjust the machine immediately incase of deviation found in products and request by quality controller for adjustment.

Prepared by: QA Coordinator  
M.Azeem

Approved by: Q.C.M  
Mr. Javeed Gafoor

Controlled Copy

**Andi Mais Co., for Medical Products**  
**Work Instruction**

**Title**  
**Production & Test Procedure for Guedel Air Way**

|         |   |             |
|---------|---|-------------|
| Doc. #  | : | WI-8-2-2-27 |
| Rev. #  | : | 7           |
| Date. # | : | 29.08.2020  |
| Page. # | : | 3 / 5       |

- 5.1.7 Quality controller is responsible to check the products every 2 hours & record the results in process control form.
- 5.1.8 Machine operator is responsible to arrange good products in clean polyethylene bag.
- 5.1.9 Production controller are responsible to weight the products and make inspection tag of products then the QC must be stamp on inspection tag to be confirm acceptable before transfer to store.

**5.2 Guedel Airway specification**

The Guedel air way dimension control / specification:

| Size | Length in (mm) | I.D. in (mm) | Colour |
|------|----------------|--------------|--------|
| 000  | 35             | 5.5          | Pink   |
| 00   | 53             | 4.5          | Blue   |
| 0    | 60             | 4.5          | Black  |
| 1    | 75             | 5.0          | White  |
| 2    | 85             | 5.5          | Green  |
| 3    | 95             | 5.5          | Orange |
| 4    | 100            | 5.5          | Red    |
| 5    | 120            | 9.5          | Yellow |

**6. Assembly procedure (By Manual)**

- 6.1 Production controller is responsible to issue parts of Guedel airway as per required size and issue production order for assembly.
- 6.2 Production controller is responsible to make sure that received parts are same as per issue slip required size of Guedel airway.
- 6.3 Assembly personnel are responsible to clean table and parts bowl with 70% IPA before assembly.
- 6.4 Guedel airway will be assembled by tight insertion of bit block in airway body without using glue.
- 6.5 Insert bite block in Guedel airway way body by tight insertion and make sure that bite block is properly inserted in airway body and not come out.

Prepared by: QA Coordinator  
M.Azeem

Approved by: Q.C.M  
Mr. Javeed Gafoor

**Controlled Copy**

**Indi Mais Co., for Medical Products**  
**Work Instruction**

**Title**  
**Production & Test Procedure for Guedel Air Way**

|        |   |             |
|--------|---|-------------|
| Doc. # | : | WI-8-2-2-27 |
| Rev. # | : | 7           |
| Date.# | : | 29.08.2020  |
| Page.# | : | 4 / 5       |

- 6.6 Quality controller is responsible to check the products during assembly as per in process control requirements of Guedel airway to make sure that products properly assembled and meeting required quality inspection criteria.
- 6.7 During assembly quality controller must be check he following visual inspection points.
- 6.7.1 The Guedel air way must be free from any irregularities defect, flashing, sink mark, gating, short shot, black spot, and blockage because of flashing. etc.
- 6.7.2 The bite block must be free from flashing, irregularities, long gating, crack and short shot, etc.
- 6.7.3 The length should be accurate as per our required specification and size.
- 6.7.4 The bite block fitting with Guedel airway should be tight insert, no loose no pull out.
- 6.7.5 Must be clean and free from dust particles and black spots.
- 6.8 Production controller is responsible to arrange the products and make inspection tag for products, then Quality controller must be stamp on inspection tag as "APPROVED".
- 6.9 After complete assembly then approved products will be transfer for packaging.
7. **Packaging and Labeling:**
- 7.1 Production controller is responsible to issue packaging film and paper as per approved packaging film and paper specification and issue production order for packaging.
- 7.2 Machine operator is responsible to set the required packaging machine die and set approved standard parameters for Guedel airway.
- 7.3 Machine operator is responsible to set the approved clichés and make sure that clichés information are correct and meeting the requirements of packaging order.
- 7.4 Machine operator is responsible to give initial sample to quality controller for inspection and get approval from quality controller before packaging of products.
- 7.5 Quality controller is responsible to the check, Product name and code no., reference no., size, lot no., type, temperature logo, sterilization logos , production date, expiry date, package forming, printing, sealing and cutting etc and quality must be stamp on pouch as approved.
- 7.6 After packaging then it will pack in specified small boxes, and arrange in big carton.
- 7.7 The big carton must be with label and detail information.

Prepared by: QA Coordinator  
M.Azeem

Approved by: Q.C.M  
Mr. Javed Gafoor

## Medi Mais Co., for Medical Products Work Instruction

### Production & Test Procedure for Guedel Air Way

|         |   |             |
|---------|---|-------------|
| Doc. #  | : | WI-8-2-2-27 |
| Rev. #  | : | 7           |
| Date. # | : | 29.08.2020  |
| Page. # | : | 5 / 5       |

#### 8. Sterilization:

8.1 The finished products will be sterile as per work instruction WI-7-5-1-1.

#### 8. Final Inspection and Product Release:

- 8.1 The Quality controller must be check, the unloaded products after sterilization.
- 8.2 The quality assurance manager or quality control manager have to review the completed cycle parameters and inspection reports both documents have to be sign.
- 8.3 After confirmation all the tests from the Lab, to release the sterilized products, the document once again have to review the quality control manager then sign.
- 8.4 The production department responsible to record the quantity of products before transferring to the finished goods ware-house and then it could be released to the customer.

#### 9. Reference:

BS EN 12181

#### 10. Sampling:

- 10.1 Inspection & test Report for Raw material
- 10.2 In-process control for Packaging Machine and final checking : F-823-T1
- 10.3 Final inspection and release record : F-823-i4  
: F-8-2-4

#### 11. Modification / Documentation change History:

| Date       | Changes                               |
|------------|---------------------------------------|
| 23.11.2018 | Item 4.3, 5.1, 5.1.1, 5.2, 8, updated |
| 29.08.2020 | Procedure 5, 6, & 7 updated           |

Prepared by: QA Coordinator  
M Azam

Approved by: Q.C.M  
Mr. Javeed Gafoor